

# DATASHEET

## QuickFix Detectable

### 1. Identification of the substances/preparation and of the company/undertaking

Date issued	09.05.2011
Revision	09.05.2014
Product name	Quick Fix Detectable
Index no.	Version 4.0
Use of the substance/preparation	Blue coloured and metaldetectable plaster for minor wounds

#### Manufacturer

Company name	Plum A/S
Postal address	Frederik Plums Vej 2
Postcode and City	5610 Assens
Country	Danmark
Tel. / Fax	+45 6471 2112 / +45 6471 2125
E-mail	<a href="mailto:info@plum.dk">info@plum.dk</a>
Website	<a href="http://www.plum.dk">http://www.plum.dk</a>

### 2. Backing

Composition	CD elastic fabric (70% Viscose, 30% polyamide)
Colour	Blue
Weight (g/m <sup>2</sup> )	90-100
MD Tensile strenght (N/cm)	>75

### 3. Adhesive

Composition	Hypoallergenic acrylic polymer
Weight (g/m <sup>2</sup> )	65-70

### 4. Coated backing

Weight (g/m <sup>2</sup> )	155-170
Adhesion to steel by peeling at 180° (N/cm)	>1,2
M.V.T.R. (g/m <sup>2</sup> /24h/37°C)	>500
Water repellency / resistance to wettability (%)	> 90

### 5. Absorbent pad

Total weight (g/m <sup>2</sup> )	225 - 275
Absorbency (g/100 cm <sup>2</sup> )	5-10
Detectable layer	Aluminium foil facing the backing
Absorbent layer	Viscose, polypropylene/polyethylene nonwoven
Colour	White
Antiadherent layer	Microperforated PE film

## 6. Liner

### Monosiliconated paper

Weight (g/m <sup>2</sup> )	57-63
Colour	White

## 7. Primary packing

	Individual pouch: Cold sealed paper + clear film
Paper Composition	Cold sealed paper
Paper weight (g/m <sup>2</sup> )	40-50
Film composition	Polyethylene terephthalate / polyethylene
Film weight (g/m <sup>2</sup> )	52-67

## 8. Secondary Packaging and outer packing

Refill Cardboard Wrap for 45 plaster units	Cardboard weight 9 gram
Inner Carton for 6 Refill Units	Cardboard weight 61 gram
Outer Carton for 8 Inner Cartons	Cardboard weight 193 gram

## 9. Handling, storage and instructions for use

Further information	To be used before expiry date. Storage temperature: ambient For single-use.
	QuickFix Detectable is a medical devices according to <i>Council Directive 93/42/EEC of 14 June 1993 concerning medical devices</i> .

## 10. Document changes

10-08-2012 /KIS	Now water repellent backing and primary packing paper
09-05-2014/KIS	Clear film added to primary packing